



Discovery* MR750

3.0T High-resolution Carotid Plaque MR

With 3.0T, 6-channel Carotid Coil and MR PlaqueView Analysis

By Kevin DeMarco, MD, Associate Professor, Director of Magnetic Resonance Imaging
Department of Radiology, Michigan State University, US

Carotid atherosclerosis may be responsible for as much as one-third¹ of all strokes. To better gauge risk of stroke from carotid atherosclerosis, clinicians are beginning to look beyond traditional risk factors such as carotid artery narrowing (stenosis) to identify high-risk features of the atherosclerotic plaque itself.^{1, 2, 3} MR plaque characterization can play a key role identifying vulnerable plaques in these patients, especially those with transient ischemic attack or in the sub-acute stage.

Patient history

A 72-year-old with asymptomatic right carotid stenosis. The patient has a history of hyperlipidemia, hypertension, and tobacco use, and has a BMI of 30.2. The patient's hyperlipidemia is well controlled (LDL-C = 88 mg/dL) with statin use since 2005. Despite meeting current AHA cholesterol guidelines, the patient has a moderately sized deep necrotic core with hemorrhage involving the right carotid bifurcation. MR carotid plaque imaging scans were performed in 2010 and in 2011 as a follow up.




Figure 1. MRI-PlaqueView delineates plaque components such as calcium (white), necrotic core (yellow), and hemorrhage (blue). (A) T1w (B) T2w (C) TOF.




Figure 2. Hemorrhage area (blue) expanded from 2010 (A) to 2011 (B).



Figure 3. Necrotic core area increased.



Kevin DeMarco, MD

is an Associate Professor and Director of Magnetic Resonance Imaging, Department of Radiology, Michigan State University, East Lansing, MI.

MR Parameters

Sequence:	Non-contrast MRA 3D TOF	Black-blood T1w 2D DIR FSE	Black-blood T2w 2D FSE
Scan time:	4:46	7:00	6:18
TR:	23	750 (1RR)	1500 (2RR)
TE:	3.9	6	48
Flip:	20	N/A	N/A
BW:	31.2	41.7	41.7
ETL:	N/A	12	16
FOV:	16	16	16
Slice thickness:	2 mm/1 mm ol	2	2
Frequency:	288	256	256
Phase:	256	224	224
Freq dir:	L/R	L/R	L/R
Phase FOV:	1	1	1
NEX:	1	2	1
# RR:	N/A	1	2

MR technique

3D TOF MRA, T1w Black-blood DIR FSE, and T2w Black-blood DIR FSE were performed using a 6-channel carotid coil. 3D TOF was used for assessing degree of stenosis, while T1w and T2w imaging were used to differentiate plaque components.

MR findings

The total volume of necrotic core was estimated at 30.1% volume in 2010 with 9.3% volume of hemorrhage. In 2011, his moderately sized hemorrhage now extended from deep in the plaque to be juxtaluminal with a new rupture of his fibrous cap. The lumen of the right ICA went from a 50% diameter stenosis by NASCET criteria with irregular narrowing to a 58% stenosis and a new ulcer in 2011. There was a much smaller hemorrhagic necrotic core involving the left ICA (not shown) with no change in the percentage of stenosis (21% in 2010 and 24% in 2011).



DIGITAL DIVE

Click to view remaining parameters on-line.




Figure 4.

Discussion

This case illustrates the potential for carotid plaque imaging to give the clinician additional information to stratify risk and help guide medical therapy. Despite very good therapeutic control of the patient’s hyperlipidemia, there is evidence of ongoing carotid plaque disease with hemorrhagic necrotic core and an increase in the right carotid stenosis, as well as the development of a new ulceration. The additional information provided by 3.0T MR carotid plaque imaging could help identify patients with ongoing plaque pathology, despite adequate therapy by current standards that may benefit from new medical treatments. **S**

3.0T Carotid Coil

A 6-channel array (3-channels per side) dedicated for high-resolution carotid imaging. It features a concentric design for deep tissue penetration and excellent SNR. ASSET is supported for accelerated imaging.

MRI-PlaqueView from VP Diagnostics

Provides visualization and quantitative analysis of multi-contrast carotid plaque MR data. Helps delineate vessel wall, calcifications and soft plaque.

References:

1. Altaf N, Daniels L, Morgan PS, et al. Detection of intraplaque hemorrhage by magnetic resonance imaging in symptomatic patients with mild to moderate carotid stenosis predicts recurrent neurological events. *J Vasc Surg.* 2008;47:337–42.
2. Takaya N, Yuan C, Chu B. Association between carotid plaque characteristics and subsequent ischemic cerebrovascular events: A prospective assessment with MRI—initial results; *Stroke* 2006; 37:818-23.
3. Singh N, Moody A R, Gladstone B J. Moderate Carotid Artery Stenosis: MR Imaging-Depicted Intraplaque Hemorrhage Predicts Risk of Cerebrovascular Ischemic Events in Asymptomatic Men. *Radiology* 2009; 252:502-8.

Kevin DeMarco, MD, is an Associate Professor and Director of Magnetic Resonance Imaging, Department of Radiology, Michigan State University, East Lansing, MI. Dr. DeMarco received his medical degree from St. Louis University Medical School and completed his internship and residency at the Naval Hospital Oakland, Oakland, CA and a neuroradiology fellowship at the University of California, San Francisco, CA. At Michigan State University, Dr. DeMarco directs a non-invasive in vivo 3.0 T MR carotid plaque research program to evaluate patients referred to specialists for known carotid stenosis. He has authored numerous publications on carotid imaging and owns two patents; one on Methods and Apparatus for Defining a Three-Dimensional Imaging Section and the other on Apparatus and Method for Detecting and Classifying Atherosclerotic Plaque Hemorrhage.

Michigan State University (MSU) is one of the top research universities in the world—on one of the biggest, greenest campuses in the nation. Home to nationally ranked and recognized academic, residential college, and service-learning programs, MSU is a diverse community of dedicated students and scholars, athletes and artists, scientists and leaders.

The MSU Department of Radiology advances knowledge and transforms lives by providing expert patient-centered medical care, delivering exceptional pre-professional and professional medical education, and expanding the understanding of wellness and disease through research and community service.